





New Voyage




ENGLISH • FRANÇAIS • DEUTSCH • ESPAÑOL • PORTUGUÊS

A

ENGLISH pag. 5

FRANÇAIS pag. 11

DEUTSCH pag. 18

ESPAÑOL pag. 25

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New Voyage

 Please read the instructions and the important notes (1.3) carefully before using the device.

1 - IMPORTANT INDICATIONS

1.1 - APPLICATION

NEW VOYAGE is an aerosol therapy device for administering drugs in aerosol form.

NEW VOYAGE incorporates a series of innovative solutions that allow rapid and effective topical treatment of the lower airways.

The MEFAR2000 kit (nebuliser + mouthpiece) or paediatric/adult mask are effective in treating disorders of the lower airways (including bronchial asthma, chronic bronchitis, pulmonary emphysema and cystic fibrosis).

 The administration of drugs requires a medical prescription specifying the type of drug, the dosage and the length of treatment.

This aerosol therapy device and its accessories have been designed and made in Italy in accordance with EU Directive 93/42/EC regarding medical devices (and subsequent updates) and bear the CE mark.

This product complies with the restriction on the use of certain hazardous substances in electric and electronic devices (European Directive RoHS 2011/65/EU).

Device suitable for continuous operation.

1.2 - SYMBOLS USED

Class II device

Device with type BF applied part

Important - read the instructions carefully

Alternating current

Device ON "I"

Device OFF "O"

The device must be dismantled at the end of its lifecycle and disposed of in special bins. Do not dispose of in ordinary town waste bins (EU Directive 2002/96/EC)

IP21 Device protected against access to hazardous parts with a finger and against vertically falling water drops.

Keep the device dry

Manufacturer

Does not contain phthalates

No traces of natural latex

1.3 - IMPORTANT NOTES

- Read the instructions carefully and keep them in a safe place for future reference.
- Do not use the device or accessories for anything other than the intended purpose.
- Children must only be allowed to use the device under the supervision of an adult who has read and understands the instructions.
- Some components are small enough to swallow and may lead to suffocation.
- The power cable and connecting pipe are long enough to pose a strangulation hazard.
- Read the list of components in the table in section 7 to see if there are any components that have caused an allergic reaction in the past.
- Do not use essential oils (e.g. with menthol or eucalyptus) as they are not compatible with the material of which the nebuliser is made.
- During the session, the device must stand on a firm flat surface.
- Never cover the air vents during operation.
- Keep well away from sources of heat.
- Do not touch the plug or the electrical components wet or with damp hands.
- IP21 Device protected against access to hazardous parts with a finger and against vertically falling water drops.
- Do not pour water or other liquids onto the device. Prevent water and other liquids from getting in through the air vents.
- Do not use in the bath or shower.
- If the device should fall into liquid, unplug before removing it. Do not use the device. Send it to an approved service centre for an overhaul.
- Beware of compressed air coming out of the air outlet as this may cause a hazard; it must only be used to supply the nebuliser.
- The device is not suitable for use in the presence of anaesthetic mixtures with air, oxygen or nitrogen protoxide.
- Performance of the device is only guaranteed if all the accessories used are original ones. Only original spare parts must be used.
- The manufacturer declares that the device has a useful life of at least 5 years.
- Insert the power cord fully into the unit.
- Remove the power cord from the device before placing it into the box.

1.4 - DESCRIPTION OF THE DEVICE

The NEW VOYAGE aerosol therapy device is made up of a compressor driven by an electric motor in the basic unit, a MEFAR2000 kit (nebuliser + mouthpiece) and two masks (which can be mounted on the nebulizer instead of the mouthpiece).

1.5 - PATIENT GROUPS

The NEW VOYAGE aerosol therapy device provides highly effective inhalation therapy for patients of any age, from young children to the elderly.

1.6 - CONTENTS

The package contains (Fig. A):

the basic unit (1), a power cable (2), an accessory storage compartment (3), an ON/OFF switch (4), a filter slot and filter (5), a power socket (6), an air intake (7), a cover (8) and a nebuliser support (9).

It also contains:

- a Mefar2000 kit containing the following:
 - a nebuliser (upper section (10), atomiser (11), lower section (12))
 - a mouthpiece (13)
- two masks: a paediatric mask (14) and an adult mask (15)
- a connecting tube (16)

2 - PREPARATION AND USE

2.1 - PREPARING THE BASIC UNIT

Before use, make sure the voltage indicated on the rating plate underneath the device corresponds to the mains voltage and the wiring system complies with the current regulations.

Do not use if the cable is damaged.

Check that the casing is integral and there is no sign of dents or damage.

To access the storage compartment, raise the front of the cover to release the snap lock.

Make sure the switch (4) is in position 0 (OFF). Plug the power cable into the socket in the apparatus, insert the electric plug in the mains socket, move the switch in position I (ON) and check that air comes out of the air intake (7).

Then move the switch in position 0 (OFF).

2.2 - PREPARING AND USING THE MEFAR2000 KIT

Please read section 3.2 before using the device.

Proceed as follows to use the Mefar2000 kit:

- Make sure the atomiser (11) is correctly inserted in the central sleeve in the lower section (12) of the nebuliser.
- The lower section (12) has a graduated scale that shows the approximate amount of liquid inside. The purpose is merely to give an approximate indication of the contents. Use a graduated syringe to get an accurate dose. Pour the liquid into the tank (Fig. B), keeping to the dosage specified by your doctor.

Readings must be taken on the graduated scale with the atomiser mounted on the lower section.

- Screw the upper section (10) of the nebuliser onto the base.
- Fit the mouthpiece (13) onto the upper section (10) of the nebuliser (Fig. C).
- Move the switch (4) in position I (ON).
- During treatment, hold the nebuliser upright in your hand to enable it to work efficiently.
- If the above operations have been performed correctly, a spray will come out of the mouthpiece.

The mouthpiece must be held in the mouth, between the teeth with the lips closed, even when you breathe out. Breathe in through the mouth and out through the nose.

During use, do not obstruct the air outlet hole in the base.

When you have used up all the liquid or wish to interrupt treatment, press (0) (OFF) on the switch and place the nebuliser back on its support (9). A small amount of liquid may remain in the nebuliser at the end of the session but this has no importance.

If a small amount of moisture remains inside the air tube at the end of the session, leave the device switched on and allow it to dry by ventilation.

The presence of liquid inside the tube is the result of water vapour forming under certain environmental conditions.

2.3 - USING THE MASK

Please read section 3.3 before using the device.

To use the mask, fit it onto the upper section of the nebuliser (Fig. C).

The mask must be placed gently over the mouth and nose. The treatment consists of inhaling the atomised liquid slowly and deeply. The exhaled air is released through vents in the sides of the mask.

During use, do not obstruct the air vents at the sides of the mask.

Atomisation is very rapid and the session only lasts a few minutes.

mask are made. They may interact and cause damage (the materials are shown in section 7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS).

3 - CLEANING, DISINFECTION AND STERILISATION

Before proceeding, unplug the device and disconnect the nebuliser and tube.

If it is necessary to clean, disinfect or sterilise this medical device using a method other than that given in this section, contact the manufacturer first.

3.1 - CLEANING, DISINFECTING AND STERILISING THE BASIC UNIT

3.1.1 - CLEANING THE BASIC UNIT

The basic unit (1) must not be cleaned in water as it could get inside and damage the electrical components, which would cause a safety hazard. Clean the basic unit (1) with a soft damp cloth, with a small amount of mild (non-abrasive) detergent, if necessary.

3.1.2 - DISINFECTING THE BASIC UNIT

The basic unit can be disinfected using a soft sponge or cloth and some cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

Make sure the solution is compatible with the material of which the basic unit is made. Refer to section 7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS.

3.1.3 - STERILISING THE BASIC UNIT

The basic unit (1) must not be sterilized.

3.2 - CLEANING, DISINFECTING AND STERILIZING THE MEFAR2000 KIT

To prevent the risk of microbial contamination, the nebuliser and mouthpiece must be cleaned, disinfected, and if necessary sterilized, after each treatment and before use, following the instructions carefully.

Before, cleaning, disinfecting or sterilizing any of the components, disassemble them as follows:

- Remove the tube (16)
- Remove the mouthpiece (13) from the nebuliser
- Unscrew the upper section (10)
- Remove the atomiser (11)

3.2.1 - CLEANING THE MEFAR2000 KIT

The disassembled nebuliser (components 10, 11, 12) and mouthpiece (13) can be washed in water at an approximate temperature of 40 °C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water.

When all the parts have been cleaned, dry them with a soft cloth and place in the storage compartment.

3.2.2 - DISINFECTING THE MEFAR2000 KIT

After washing, the disassembled nebuliser (components 10, 11, 12) and mouthpiece (13) can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite). Follow the manufacturer's instructions.

The components of the nebuliser and the mouthpiece can be disinfected by boiling in a saucepan of water for a maximum of 10 minutes. Make sure the components are not in direct contact with the bottom of the saucepan.

Boiling may alter the colour or appearance of the components, but this will not affect operation of the device.

3.2.3 - STERILIZING THE MEFAR2000 KIT

All components of the nebuliser (10, 11 and 12) and the mouthpiece (13) can be sterilized in an autoclave at 121 °C for at least 15 minutes (but no more than 30 minutes), for up to 20 times.

Sterilisation is recommended before each session when the device is used in a hospital environment or by more than one patient.

3.3 - CLEANING, DISINFECTING AND STERILIZING THE MASKS

To prevent the risk of microbial contamination, the masks (14, 15) can be cleaned and disinfected, and if necessary sterilised, after each treatment and before use, following the instructions carefully.

3.3.1 - CLEANING THE MASKS

The masks can be washed in water at an approximate temperature of 40 °C with the addition of detergent (refer to the manufacturer's instructions), then rinsed in running water.

After cleaning, dry with a soft cloth and place in the storage compartment.

3.3.2 - DISINFECTING THE MASKS

After washing, the masks can be immersed in a cold disinfectant solution (e.g. a solution containing up to 2% sodium hypochlorite), following the manufacturer's instructions.

⚠ The masks (14, 15) must not be boiled.

3.3.3 - STERILIZING THE MASKS

⚠ The masks must not be sterilised hot or in an autoclave.

3.4 - CLEANING, DISINFECTING AND STERILIZING THE TUBE

The tube (16) needs to be cleaned after each session. Detach it from the nebuliser and leave the device on until all signs of moisture have been eliminated. If the tube is particularly dirty, it needs to be replaced.

⚠ The tube (16) must not be sterilized hot in an autoclave or boiled.

4 - MAINTENANCE

The NEW VOYAGE aerosol therapy device does not require any particular maintenance. The compressor does not need to be lubricated as the piston is the dry-operating type.

A special external filter keeps out dust and moisture, and should be replaced periodically.

⚠ Always unplug from the mains before carrying out the following operations.

4.1 - REPLACING THE AIR FILTER (Fig. D)

The air filter needs to be replaced once a year if the device is used frequently.

To replace the air filter (5), remove it using a screwdriver and insert a new one.

The filter cannot be washed or regenerated. It must always be replaced with a new one.

4.2 - MAINTENANCE OF THE ACCESSORIES

When cleaning, disinfecting or sterilising the components, you can check them visually and determine when they need to be replaced. Repeated sterilisation may affect the materials of which the accessories are made, so the components should be inspected regularly. The accessories must be replaced when there are clear signs of deterioration.

5 - PROBLEMS, POSSIBLE CAUSES AND SOLUTIONS

TROUBLE-SHOOTING CHART

Problems	Possible causes	Solutions
The device does not work	<ul style="list-style-type: none"> - The plug is not inserted in the socket. - The power connector is inserted wrongly. - The switch (4) is in position 0 (OFF). - The thermostat has cut in. 	<ul style="list-style-type: none"> - Insert the plug in the socket and move the switch in position I (ON). - Insert it correctly. - Move the switch in position I (ON). - Clear the air vents.
The nebuliser does not generate a spray	<ul style="list-style-type: none"> - The device is switched off. - There is no liquid in the nebuliser. - Too much liquid has been poured into the nebuliser. - The nebuliser has not been assembled correctly. - The air outlet in the base (12) is blocked. - The air tube (16) is not connected securely. 	<ul style="list-style-type: none"> - Switch on the device. - Pour liquid into the nebuliser. - Reduce the quantity of liquid in the nebuliser to below the maximum level. - Check that the atomiser (11) is properly inserted in the base of the nebuliser (12). Screw the upper section (10) onto the base and check that air comes out of the little hole at the top of the base. - Unblock the hole. - Connect the tube (16) securely to the air intake on the device (7) and the one on the nebuliser.

If you have failed to solve the problem after performing the above checks and the device still does not work, have it checked by an authorized service centre.

6 - TECHNICAL DATA


The operating characteristics vary with the version of aerosol therapy device.

As you can see in the table below, each version has its own supply voltage, which is specified on the sticker under the basic unit.

Supply voltage	NEW VOYAGE-1 CE 0051 220/240 V ~ 50 Hz	NEW VOYAGE-2 CE 0051 220 V ~ 60 Hz	NEW VOYAGE-3 240 V ~ 50/60 Hz	NEW VOYAGE-4 115 V ~ 60 Hz	NEW VOYAGE-5 110 V ~ 50/60 Hz
Maximum air pressure	2,5 bars	2,5 bars	2,5 bars	2,5 bars	2,5 bars
Maximum air flow	10 l/min	10 l/min	10 l/min	10 l/min	10 l/min
Dynamic air flow	7 l/min	7 l/min	7 l/min	7 l/min	7 l/min
Operation	Continuous				
Maximum input	See rating plate on the device				
Safety fuse (inside)	See rating plate on the device				
Dimensions of basic unit	207x285x107 (h) mm				
Weight (excluding power cable)	1.5 kg				
Noise level at 1 m	57.3 dBA (according to UNI EN 13544 -1)				
Minimum nebuliser filling volume	1 ml				
Maximum nebuliser filling volume	8 ml				
Aerosol output with nebuliser (*)	0.75 ml (**)				
Minimum aerosol output rate with nebuliser (**)	0.30 ml/min				
MMAD, Graph 1 (*) (***)	2.1 µm (**)				

CE devices bear the CE mark on the sticker under the basic unit

GRAPH 1

The MMAD values shown are measured using a MACH 2 API AEROSIZER[®].

The y-axis of the graphs represents the percentage of volume of atomised solution with sizes lower than the corresponding dimensional value expressed in µm on the x-axis.

(*) For a quantity of 2 ml.

(**) The values shown refer to the use of physiological solution (0,9% NaCl) and may vary with the drug used. The values given do not apply to drugs in suspension or high-viscosity drugs. The manufacturer of the drug can provide the necessary details.

(***) 50% by volume of the atomised particles is less than the MMAD (Median Mass Aerodynamic Diameter). This value is indicative of the nebuliser's ability to produce particles suitable for a specific therapeutic treatment.

7 - DISPOSING OF THE DEVICE AND ITS COMPONENTS

At the end of their useful life, the device and its components must be disposed of in accordance with the current regulations and environmental laws. If no legal obligations exist, the various materials must be sorted and disposed of separately.

The table below shows the materials of which the various components are made.

None of the materials used contains phthalates or traces of natural Latex.

COMPONENT	MATERIAL
Basic unit (1)	Waste from Electric and Electronic Equipment (WEEE)
External casing of the basic unit	ABS
Nebuliser (10, 12)	Polycarbonate
Atomiser (11)	Nylon
Mouthpiece (13)	Polycarbonate
Masks (14, 15)	Phthalate-free polyvinyl chloride (PVC)
Connecting tube (16)	Phthalate-free polyvinyl chloride (PVC)
Box	Cardboard
Instructions	Paper
Accessory bag	Low-density polyethylene (LDPE)

8 - ELECTROMAGNETIC COMPATIBILITY

This device complies with the requirements of CEI EN 60601-1-2:2015 on the electromagnetic compatibility of medical devices. Compliance with electromagnetic compatibility standards does not guarantee total immunity of the product. Some devices (mobile phones, pagers, etc.) may interrupt operation if used near medical devices. None of the materials used contains phthalates or traces of natural Latex.

RECOMMENDED DISTANCES BETWEEN RADIO-COMMUNICATION DEVICES

The NEW VOYAGE aerosol therapy device is designed to operate in an electromagnetic environment where radiated RF disturbance is under control. The customer or user of the NEW VOYAGE aerosol therapy device can help prevent electromagnetic interference by leaving a minimum distance between mobile and portable RF communication devices (transmitters) and the NEW VOYAGE aerosol therapy device. The table below shows the recommended distances for the NEW VOYAGE aerosol therapy device based on the maximum output of the radio communication devices.

Nominal maximum transmitter output [W]	Separation distance [m] at the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \times \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \times \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \times \sqrt{P}$
0.01	0,12	0,12	0,23
0.1	0,38	0,38	0,73
1	1,2	1,2	2,33
10	3,8	3,8	7,3
100	12	12	23

The recommended separation distance d in metres [m] for transmitters with a nominal maximum transmitter output not shown in the table can be calculated using the equation applied to the transmitter frequency, where P is the nominal maximum transmitted output in Watt [W] declared by the manufacturer.

Note 1. For 80 MHz and 800 MHz, apply the distance for the higher frequency interval.
Note 2. These guidelines may not apply in all situations. Electromagnetic propagation is influenced by the absorption and reflection of structures, objects and persons.

9 - OPERATING AND STORAGE CONDITIONS

9.1 - OPERATING CONDITIONS

- Temperature range +5 °C to +40°C
- Relative humidity 15% to 93%
- Atmospheric pressure 700 hPa to 1060 hPa

9.2 - STORAGE CONDITIONS

- Temperature range -25 °C to +50 °C
- Relative humidity less than 93%
- Atmospheric pressure 700 hPa to 1060 hPa

10 - SPARE PARTS

Only original spare parts must be used.

11 - WARRANTY CONDITIONS

The manufacturer is liable for the safety, reliability and performance of the device, provided that it is used in accordance with the instructions and for the intended purpose only, and any repairs are carried out by the manufacturer or an approved service centre. The manufacturer declines all liability and the warranty will be invalidated if the user fails to follow the instructions and the important notes contained herein.

12 - TECHNICAL UPDATES

Air Liquide Medical Systems S.r.l. periodically reviews and modifies all its medical devices in order to improve their performance, safety and reliability. The instruction booklets are updated accordingly to include any new or changed features. If the booklet accompanying the device gets damaged or mislaid, a replacement copy can be obtained from the manufacturer by quoting the data shown on the rating plate.

13 - COPYRIGHT

The information contained in this booklet must not be used for any purpose other than the reason for which it is provided. This booklet is the property of Air Liquide Medical Systems S.r.l. and no part of it may be reproduced without the manufacturer's written permission. All rights reserved.

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